

Plethysmography is used to measure respiratory functions in both human and lab animal. Respiratory waveforms, air analysis, temperature and electrophysiology are recorded into data sets. Currently, we use LabChart to record data. We will pick up data set manually and transfer the data to excel sheet and MATLAB for further analysis. This process is time consuming and artificial.

To get more consistent results, I use "Plethysmography.py" to automated analyze multiple input files (.txt) and produce figures.

To run the "Plethysmography.py", simply enter the folder of the program by type "cd <path." and the type "python ./Plethysmography.py".

Experimental conditions are required by user inputs and are saved as "MUID.csv" under subfolder "/ExpInfo".

Experimental data inputs are converted by LabChart to "MUID.txt" and "MUID_comment.txt". (MUID indicate the experimental ID of each mouse)

Each mouse data will be saved as "MUID_data.csv" and will be summarized into "MUID.csv" under subfolder "/MUID"


Figures will be then generated and saved as "Title.tif" under subfolder "/figures")

Title includes 8 different parameters we are looking at: Respiratory Rate, Tidal Volume, Tidal Volume (Body Weight normalized), Minute Ventilation, Oxygen Consumption, Carbon Dioxide, Minute Ventilation normalized to Oxygen Consumption, and Temperature.


Packages required: os, numpy, scipy, pandas, matplotlib, datetime, csv, neurokit

Examples:

Input:



Output:



```
main()
Start a new project? (y/n): y

ExpCon, Treatment, GenoGroups, Dose = NewProj()

Experimental Conditions (hypoxia or hypercapnia): hypoxia
CNO? (y/n): y
Dose (mg/kg): 1
List all the genotypes in this project (separate by ;): SST_FLPo x F_hM4De;WT
```

```
path, files = LoadFiles()
```

```
The path of file:/Users/huananishi/Documents/Graduate school/Second year /Python/project  
File (MUIDs) names (separate by ;):M27699;M28041;M28042;M28043;M28045;M28641
```

```
file2, Stat = ExpInfo(files,path)
```


```
M27699  
Sex (m/f): m  
Genotype: SST_FLPo x F_hM4De  
Group (Exp or Con): Exp  
Body Weight (g):25.4  
Beginning Temperature (C): 39.2  
Middle Temperature(C): 36.4  
End Temperature (C): 36.1  
Post 30 min Temperature (C): 37.2  
Room Temperature(F): 69.5  
Barometric Pressure (inHg): 29.58  
Experiment Date (mm/dd/yy): 11/18/17  
Birth Date (mm/dd/yy): 9/22/17  
M28041  
Sex (m/f): f  
Genotype: WT  
Group (Exp or Con): Con  
Body Weight (g):21  
Beginning Temperature (C): 37.7  
Middle Temperature(C): 37  
End Temperature (C): 36.9  
Post 30 min Temperature (C): 37.7  
Room Temperature(F): 67  
Barometric Pressure (inHg): 29.7  
Experiment Date (mm/dd/yy): 1/15/18  
Birth Date (mm/dd/yy): 10/25/17
```

...

```
DAQ(files,path)
```

```
***** M27699 done*****  
***** M28041 done*****  
***** M28042 done*****  
***** M28043 done*****  
***** M28045 done*****  
***** M28641 done*****
```

```
Report(path,Stat,Treatment,ExpCon,CompiledData,CompiledTemp,GenoGroups)
```



...